

## BLOOD PRESSURE MONITOR TM-2655/TM-2655P WITH PRINTER

Fast. Easy. Accurate.



The TM-2655 Blood Pressure Monitor for health promotion and screening builds customer loyalty and store traffic. It is used in retail, medical, or occupational health settings.

The TM-2655 is a fully automatic blood pressure monitor that ensures proper measurement through proper positioning. Unlike other commercial blood pressure monitors, the TM-2655 allows the patient more flexibility in positioning when taking a measurement.

It is uniquely designed for patients to use either the left or right arm for measurement.

More importantly, the patient does not have to twist their body to insert the arm in the cuff. This feature helps people become more comfortable when taking a measurement. It is especially convenient for pregnant women, patients with back problems and patients in wheelchairs.

A large display shows the measurement results. Systolic, diastolic and pulse rate are displayed simultaneously.

This monitor is available with or without a printer. Printouts can be customized to include your store's address or advertisement.

	Feb, 11 2010 10:52
<b>SYS</b>	120 mmHg
<b>DIA</b>	80 mmHg
<b>PUL</b>	72 bpm
<b>Pharmacy Name</b>	
Street Address	
City, State Zip	
Phone #	
Message/Comments	

- Compact, tabletop design
- Can use right or left arm for measurement
- One button operation
- Quick printout (printer version only)
- Large display
- Displays date and time
- Energy saving function
- Clinically validated

# TM-2655/TM-2655P Blood Pressure Monitor

You can select a TM-2655 unit with or without a printer. The TM-2655P unit with a printer can be customized before it is shipped to you so that it can include your address or advertisement on the blood pressure measurement printouts.

PART NUMBER	UPC NUMBER	DESCRIPTION
TM-2655	0-93764-60182-8	Without Printer
TM-2655P	0-93764-60062-3	With Printer

Each unit comes with a power cable, 1 extra cuff cloth, instruction manual and instruction card with plastic holder. TM-2655P comes with 2 rolls of printer paper.

## Accessories

PART NUMBER	DESCRIPTION
AX-133003442-S	TM-2655 Cuff Cloth (5 pieces/set)
AX-KO1869	TM-2655 RS-232C Cable
08:4009441	TM-2655 Instruction Sheet
07:3000036	TM-2655 Instruction Sheet Holder
PP-751	Printer Paper (pack of 5)

## Technical Specifications

<b>Dimensions</b>	15.4 (l) x 9.6 (w) x 12.8" (h) (390 x 245 x 325 mm)
<b>Weight</b>	Approximately 19.8 lbs (9.0 kgs)
<b>Display</b>	Digital, LED - Systolic, Diastolic and Pulse
<b>Time/Date</b>	24 hour format, year (up to 2091), month and day
<b>Pressurization</b>	Automatic, using micropump
<b>Deflation Type</b>	ECEV and Automatic Constant Speed Exhaust Valve (Dual) Oscillometric
<b>Accuracy</b>	Pressure: $\pm 3$ mmHg or $\pm 2\%$ , whichever is greater; Pulse: $\pm 5\%$
<b>Measurement Range</b>	Pressure: 0 mmHg to 300 mmHg; Pulse: 30 to 200 pulses/minute
<b>Power Source</b>	100/115, 220/240VAC, 50/60Hz
<b>Printer</b>	Thermal Type, 58 mm wide, Automatic Cutter
<b>Safety Device (Electrical)</b>	Quick Release when START button is pressed Quick Release when EMERGENCY Stop Button is pressed Automatic Rapid Deflation when pressure exceeds 320 mmHg Automatic Rapid Deflation / Power Off when pressure reaches 320 to 350 mmHg
<b>Safety Device (Mechanical)</b>	Releases cuff when safety lever is pressed
<b>Operating Environment</b>	50°F to 104°F (10°C to 40°C), 30% to 85% RH
<b>Storage Environment</b>	-4° to 140° F (-20° to 60° C), 10% to 95% RH

Blood pressure measurements determined with this device are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultation method within the limits prescribed by the American National Standards Institute (ANSI/AAMI SP10) for electronic sphygmomanometers. Manufactured by an ISO 9001 and ISO 13485 registered company. Engineered and made in Japan.

## Packaging and Shipping Specifications

<b>Selling Unit (Box)</b>	<b>L x W x H (inches)</b>	<b>WT (lb.)</b>
TM-2655	20 <sup>5</sup> / <sub>8</sub> x 14 <sup>1</sup> / <sub>8</sub> x 19	27.25
<b>Unit Case (10 units)</b>	<b>L x W x H (inches)</b>	<b>WT (lb.)</b>
TM-2655P	20 <sup>5</sup> / <sub>8</sub> x 14 <sup>1</sup> / <sub>8</sub> x 19	27.50



1756 Automation Parkway  
San Jose, CA 95131  
1-800-726-3364  
www.andmedical.com



Ask your A&D Sales Representative about other A&D Medical products.

- Blood Pressure Monitors
- Digital Thermometers
- Personal Scales
- Stethoscopes
- Pedometers

For questions about A&D Medical products, call us toll-free at 1-800-726-3364 or visit our web site at [www.andmedical.com](http://www.andmedical.com)