

## ONE-STEP PLUS MEMORY Blood Pressure Monitor UA-767PV

Fast. Easy. Accurate.



**A sleek, easy-to-use blood pressure monitor with one button operation plus advanced features.**

The UA-767 Plus blood pressure monitor incorporates the latest technology for accurate measurements at home. The One Step Plus Memory monitor is fully automatic and inflates to the correct pressure each time a measurement is taken.

This monitor offers the exclusive Pressure Rating Indicator™, Irregular Heartbeat Detection, Clock, Time & Date Stamp, 90 memory recall and is clinically validated for professional accuracy.

- **Pressure Rating Indicator™** - providing a snapshot of the user's blood pressure classification according to the US JNC7 guidelines.
- **Irregular Heartbeat Detection** - will alert the user to the presence of an irregular heartbeat while still providing an accurate reading.
- Provides simultaneous readout of systolic/diastolic pressure and pulse rate
- Calculates and displays average blood pressure readings
- Available in small, medium or large cuff sizes (cuffs available separately)
- Large-print trilingual instruction guide with illustrations and quick start card in English, French and Spanish
- Four week blood pressure logbook

For questions about LifeSource monitors, call us toll-free at 1-800-726-3364 and a specially trained representative will assist you, or visit us at [www.LifeSourceOnline.com](http://www.LifeSourceOnline.com).

- 90 Memory Recall
- Clinically validated for accuracy
- Clock, Time & Date Stamp
- Pressure Rating Indicator™
- Irregular Heartbeat Detection
- One button operation
- Displays Average Readings
- Large digital display
- Fast measurement
- Latex free
- Optional AC Power Adapter sold separately

# UA-767 Plus - One Step Plus Memory

To obtain an accurate reading, it is important to use the correct cuff size. A cuff that is too small will result in a higher than expected reading and a cuff that is too large will result in a lower than expected reading. The UA-767 Plus - One Step Plus Memory unit comes in three different cuff sizes to accommodate your customers' needs.

|               | PART NUMBER | UPC NUMBER      | CUFF SIZE                    |
|---------------|-------------|-----------------|------------------------------|
| UNITED STATES | UA-767PV    | 0-93764-60047-0 | Medium 9.4 - 14.2" (24-36cm) |
|               | UA-767PVL   | 0-93764-60048-7 | Large 14.2 - 17.7" (36-45cm) |
|               | UA-767PVS   | 0-93764-60056-2 | Small 6.3 - 9.4" (16 - 24cm) |
| CANADA        | UA-767PCN   | 0-93764-60048-4 | Medium 9.4 - 14.2" (24-36cm) |
|               | UA-767PCNL  | 0-93764-60050-0 | Large 14.2 - 17.7" (36-45cm) |
|               | UA-767PCNS  | 0-93764-60057-9 | Small 6.3 - 9.4" (16 - 24cm) |

Each unit comes with a trilingual instruction manual, quick start card, and four week blood pressure logbook.

## Accessories

|                          | PART NUMBER | UPC NUMBER      | CUFF SIZE                       |
|--------------------------|-------------|-----------------|---------------------------------|
| UNITED STATES AND CANADA | UA-279      | 0-93764-01154-2 | Small: 9.4 - 14.2" (24-36cm)    |
|                          | UA-280      | 0-93764-04607-0 | Medium: 9.4 - 14.2" (24-36cm)   |
|                          | UA-281      | 0-93764-01069-9 | Large: 14.2 - 17.7" (30 - 45cm) |
|                          | TB:233      | 0-93764-60235-1 | AC Power Adapter                |

## Technical Specifications

|                           |   |
|---------------------------|---|
| <b>Dimensions</b>         | 4.3" (l) x 5.6" (w) x 2.5" (h) (110 mm x 142 mm x 64 mm)  |
| <b>Weight</b>             | 10.6 oz. (300 g) without batteries  |
| <b>Display</b>            | Digital character height: 16 mm (blood pressure), 10 mm (pulse rate)<br>Pressure and pulse displayed simultaneously         |
| <b>Memory</b>             | 90 readings   |
| <b>Pressurization</b>     | Automatic, using micropump  |
| <b>Deflation</b>          | Automatic exhaust   |
| <b>Type</b>               | Oscillometric   |
| <b>Accuracy</b>           | Pressure: $\pm 3$ mmHg or $\pm 2\%$ , whichever is greater; Pulse: $\pm 5\%$  |
| <b>Pressure Range</b>     | 20 mmHg to 280 mmHg; Pulse: 40 to 200 pulses/minute   |
| <b>Power Source</b>       | 120V AC adapter (model TB:181A) or four type AA alkaline batteries (1.5 volt).<br>Batteries and power adapter not included. |
| <b>Battery Life</b>       | Approx. six months with one daily measurement   |
| <b>Operating Environ.</b> | 50° to 104° F (10° to 40° C)  |
| <b>Storage Environ.</b>   | 14° to 140° F (-10° to 60° C)   |

Blood pressure measurements determined with this device are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultation method within the limits prescribed by the American National Standards Institute (ANSI/AAMI SP10) for electronic sphygmomanometers. Manufactured by an ISO 9001 and ISO 13485 registered company. Engineered in Japan. Made in China.

## Packaging and Shipping Specifications

|                             |                           |                 |  |
|-----------------------------|---------------------------|-----------------|--|
| <b>Selling Unit (Box)</b>   | <b>L x W x H (inches)</b> | <b>WT (lb.)</b> | <b>Individual Item Cube (cubic feet)</b> |
| UA-767 Plus Models          | 5 1/8 x 3 15/16 x 7 11/16 | 1.49            | 0.09                                     |
| UA-767 PAC Models           | 5 1/8 x 3 15/16 x 7 11/16 | 1.56            | 0.09                                     |
| <b>Unit Case (10 units)</b> | <b>L x W x H (inches)</b> | <b>WT (lb.)</b> | <b>Case Cube (cubic feet)</b>            |
| UA-767 Plus Models          | 20 5/8 x 10 7/8 x 9       | 16.02           | 1.17                                     |
| UA-767 PAC Models           | 20 5/8 x 10 7/8 x 9       | 17.04           | 1.17                                     |



1756 Automation Parkway  
San Jose, CA 95131  
www.LifeSourceOnline.com

LifeSource is a trademark of A&D Medical. ©2010 A&D Medical. All rights reserved.  
Pressure Rating Indicator is a trademark of A&D Medical and is patent pending.



Ask your A&D Sales Representative about other LifeSource products.

- Blood Pressure Monitors
- Digital Thermometers
- Personal Scales
- Stethoscopes
- Pedometers

For questions about *LifeSource* monitors, call us toll-free at 1-800-726-3364 or visit our web site at [www.LifeSourceOnline.com](http://www.LifeSourceOnline.com)

LIT:LS767P 09/10