

PREMIUM AUTOMATIC UA-853AC Blood Pressure Monitor

Fast. Easy. Accurate.



A monitor with advanced features in a sleek, contemporary design.

The sleek new Premium Blood Pressure Monitor combines clinically valuable features with a contemporary new table top design. This monitor offers our new AM/PM TimeWise™ Tracking Technology, Pressure Rating Indicator™, Irregular Heartbeat Detection, three programmable alarms, 90 memory recall, time and date features.

- **AM/PM TimeWise™ Tracking Technology** - which automatically categorizes and averages readings as AM or PM to allow patient and health care provider to achieve optimal treatment and blood pressure goals.
- **Pressure Rating Indicator™** - providing a snapshot of the user's blood pressure classification according to the JNC7 guidelines.
- **Irregular Heartbeat Detection** - will alert the user to the presence of an irregular heartbeat while still providing an accurate reading.
- **Three programmable alarms** - to remind users to take a reading, or for use as a daily alarm.
- **90 memory** – records the last 90 readings.
- AC Adapter included
- Trilingual instruction guide with illustration and quick start card
- Four week blood pressure logbook included

LifeSource blood pressure monitors lead the way in technology and quality – because peoples' lives depend on them. Our commitment to patient-focused thinking has made us the leading primary manufacturer in developing clinically valuable features. The LifeSource blood pressure monitors are fast, easy, accurate and carry a full Lifetime Warranty.

Visit our website at www.LifeSourceOnline.com for more information or call our Customer Satisfaction hotline at 1-(888) 726.9966.

- AM/PM TimeWise™ Tracking Technology
- 90 Memory Recall
- Pressure Rating Indicator™
- Irregular Heartbeat Detection
- Displays Average Readings
- Three programmable alarms
- Large digital display
- AC Adapter included
- Fast measurement

Premium Automatic UA-853AC Blood Pressure Monitor

To obtain an accurate reading, it is important to use the correct cuff size. The UA-853AC Premium Automatic monitor includes either a Medium or a Large cuff. A cuff that is too large will produce a reading that is lower than the correct blood pressure; a cuff that is too small will produce a reading that is higher than the correct blood pressure.

	PART NUMBER	UPC NUMBER	CUFF SIZE
UNITED STATES	UA-853AC	0-93764-60194-1	9.4-14.2" (24-36cm)
	UA-853LAC	0-93764-60195-8	14.2-17.7" (36-45cm)
CANADA	UA-853CNAC	0-93764-60197-2	9.4-14.2" (24-36cm)
	UA-853CNLAC	0-93764-60198-9	14.2-17.7" (36-45cm)

Each unit comes with an AC Adapter, carrying pouch, trilingual instruction manual, quick start card, four week blood pressure logbook and warranty card.

Accessories

	PART NUMBER	UPC NUMBER	CUFF SIZE
UNITED STATES	UA-280	0-93764-04607-0	Medium: 9.4" - 14.2" (24-36cm)
	UA-281	0-93764-01069-9	Large: 14.2" - 17.7" (30 - 45cm)
	TB:181C or TB:233	NA	AC Power Adapter

Technical Specifications

Dimensions	5.9" (w) x 4.9" (l) x 6.1" (h) (150mm x 126mm x 156mm)
Weight	12.3 oz.. (350g) without batteries
Display	Digital, 20-mm character height Pressure/pulse displayed simultaneously
Time/Date	12 hour format (AM/PM), year (2006-2056), month, and day with automatic adjustment for leap years
Pressurization	Automatic, using micropump
Deflation	Active exhaust valve
Type	Oscillometric
Accuracy	Pressure: ± 3 mmHg or $\pm 2\%$, whichever is greater Pulse: $\pm 5\%$
Pressure Range	Pressure: 20 mmHg to 280 mmHg Pulse: 40 pulses to 180 pulses
Power Source	4 type "AA" (1.5 volt) alkaline batteries (not included) or 120 V AC Adapter
Battery Life	Approximately 4 months with one daily measurement
Operating Environment	50°F to 104°F (10°C to 40°C) Less than 85% relative humidity
Storage Environment	14°F to 140°F (-10°C to 60°C) Less than 95% relative humidity

Blood pressure measurements determined with this device are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultation method within the limits prescribed by the American National Standards Institute (ANSI/AAMI SP10) for electronic sphygmomanometers. Manufactured by an ISO 9000 and ISO 13485 registered company. Engineered in Japan. Made in China.

Packaging and Shipping Specifications

Selling Unit (Box)	W x D x H (inches)	WT (lb.)	Individual Item Cube (cubic feet)
UA-853AC/CNAC	6.3 x 6.7 x 7.7	2.1	0.19
UA-853LAC/CNLAC	6.3 x 7.9 x 7.7	2.1	0.22
Unit Case (10 units)	W x D x H (inches)	WT (lb.)	Case Cube (cubic feet)
UA-853AC/CNAC	14.2 x 32.3 x 9.3	23.3	2.47
UA-853LAC/CNLAC	16.4 x 32.3 x 8.5	23.3	2.61

Ask your A&D Sales Representative about point-of-purchase floor and counter displays available.



1756 Automation Parkway
San Jose, CA 95131
www.LifeSourceOnline.com

LIFE SOURCE



Pressure Rating Indicator is a trademark of A&D Medical and is Patent Pending. TimeWise is a trademark of A&D Medical. LifeSource is a registered trademark of A&D Medical. ©2008 A&D Medical. All rights reserved.



Ask your A&D Sales Representative about other LifeSource products.

- Blood Pressure Monitors
- Digital Thermometers
- Personal Scales
- Stethoscopes
- Pedometers

For questions about LifeSource products, call us toll-free at 1-800-726-3364 or visit our web site at www.LifeSourceOnline.com

LIT:LS853 02/08