

Fast. Easy. Accurate.

WELLNESS CONNECTED WIRELESS PRECISION SCALE UC-324THX



The UC-324THX Wireless Precision Scale steps up to your wellness needs.



The Wellness Connected Wireless™ Precision Scale helps provide you with a complete record of your weight and BMI. The UC-324THX features an easy-to-read LCD display and internal memory storage of 63 measurements. Also features one-touch activation, superior accuracy, high capacity, low profile form factor, and a two year warranty. Track your daily weight online at www.WellnessConnected.com

through any computer with an internet connection. View your daily, weekly or monthly weight measurements and compare it against your personal goals to keep you moving towards better health! No buttons to press, your activity data is automatically uploaded to any computer.

- **Wireless Connectivity** – Using patented wireless technology your blood pressure readings will be effortlessly transmitted to your online Wellness Connected service account.
- **Cloud Computing** – Your data is stored in the cloud, meaning you can access it from anywhere you need with an internet connection!
- **Precision Accuracy** – Measures weight in 0.2 lb (0.1 kg) increments allowing you to track your results with absolute precision.
- **Body Mass Index Calculation** – Automatically calculates Body Mass index (BMI) according to height.
- **High Capacity** – 330 lb (150 kg) capacity.
- **One-Touch Activation** – Fast and easy to get your results and keep on moving.
- **4 AA Batteries Included.**



- Wireless Connectivity
- User-friendly Software
- One Touch Activation
- Automatic Data Transfer
- Capacity of 330lbs (150kg)
- Precision Accuracy
- 60+ Measurement Internal Memory
- Trilingual Instructions
- Track Your Weight and BMI Online

For questions about Wellness Connected products, call us toll-free at 1-888-726-9966 and a specially trained representative will assist you, or visit us at www.WellnessConnected.com.

UC-324THX Wireless Precision Scale

PART NUMBER	UPC NUMBER	CUFF SIZE
UC-324THX	0-93764-60281-8	Wireless Activity Monitor

Contents Include:

- Wireless Precision Scale
- Trilingual Instruction Manual (English, French, Spanish)
- 4 AA Batteries

Wellness Connected Online Service



Join to get your real-time health and activity data at WellnessConnected.com

Technical Specifications

Measuring Method	Electric Resistance
Digital Display	Liquid Crystal Display
Maximum Capacity (MAX)	330 lbs / 150 kg
Minimum Display (d)	0.2 lbs / 0.1 kg
Power Source	Four Type AA Batteries (1.5 V)
Dimensions	13.8 (D) x 13.8 (W) x 1.5 (H) in 350 (D) x 350 (W) x 39 (H) mm
Weight	4.8 lbs (2.2 kg), including the batteries

Packaging and Shipping Specifications

Selling Unit (Box) UC-324THX	W x D x H (inches) 16.4 x 3.3 x 16.3	Wt (lb.) 6.64	Individual Item Cube (cubic feet) 0.51
Unit Case (5 units) UC-324THX	W x D x H (inches) 17.8 x 17.6 x 17.5	Wt (lb.) 36.0	Case Cube (cubic feet) 3.17

Ask your A&D Medical Sales Representative about the entire Wellness Connected family of products.

A&D
A&D Medical

LIFE SOURCE
A&D Medical



1756 Automation Parkway
San Jose, CA 95131
1-888-726-9966
www.andmedical.com

Wellness Connected and LifeSource are registered trademarks of A&D Medical. ©2011 A&D Medical. All rights reserved.



- 1 Tap scale
- 2 Step on
- 3 Read results on display or software

Ask your A&D Medical Sales Representative about other A&D Medical & Wellness Connected products.

- Blood Pressure Monitors
- Digital Thermometers
- Personal Scales
- Stethoscopes
- Activity Monitors

For questions about A&D Medical products, call us toll-free at 1-888-726-9966 or visit one of our websites at:
www.andmedical.com
www.WellnessConnected.com
www.LifeSourceOnline.com

LIT-UC-324THX 06/11