




## Stethoscopes



The **LifeSource** Stethoscope line includes Single Head, Dual Head and Two-in-One Tube Sprague Rappaport models in a variety of colors.

**LifeSource**<sup>™</sup> Stethoscopes offer practical and affordable acoustic performance to physicians, nurses, students and others. Latex-free tubing protects patients and medical professionals from allergens while providing a durable and reliable measurement tool.

The stethoscopes are constructed of fully adjustable stainless steel binaurals with ear tips designed to provide comfort and enhanced sound transmission. The rotating chest piece on the Dual Head and Sprague Rappaport Stethoscopes allows flexibility in distinguishing high and low frequencies. **LifeSource** Stethoscopes come with additional ear tips and diaphragm, and are available in four tube colors (number indicates color code): black (01), navy (05), purple (09) or hunter green (11).

- Three styles with four different colors to choose from
- Comes in retail package
- Affordable
- Latex free
- Lifetime Warranty

For questions about **LifeSource** products, call our toll free **LifeSource** HealthLine at 1-888-726-9966 and a specially trained representative will assist you, or visit us at [www.LifeSourceOnline.com](http://www.LifeSourceOnline.com).

- Patient-friendly, latex-free tubing
- Convenient 30" overall length, 22" tube length
- Durable chrome-plated binaurals
- Fully rotating chestpiece on Dual Head and Sprague Rappaport Stethoscopes
- Various colors available



# Stethoscopes

*LifeSource* stethoscopes come in three models to fit someone's needs. Each model comes in four select colors.

## Single Head Stethoscope

- Lifetime warranty
- Latex-free tubing
- 30" overall length
- 22" Y-tubing
- Chrome-plated binaurals
- Aluminum chestpiece
- Extra ear tips and diaphragm included



## Dual Head Stethoscope

- Lifetime warranty
- Latex-free tubing
- 30" overall length
- 22" Y-tubing
- Chrome-plated binaurals
- Fully rotating aluminum chestpiece
- Extra ear tips and diaphragm included



## Sprague Rappaport Stethoscope

- Lifetime warranty
- Latex-free tubing
- Advanced two-in-one design
- 30" overall length
- 22" Y-tubing
- Chrome-plated binaurals and fully rotating chestpiece
- 2 extra ear tips (mushroom and bell style)
- 3 bell rings (adult, medium and infant)
- Extra diaphragms included



COLOR	SINGLE HEAD	UPC NUMBER	DUAL HEAD	UPC NUMBER	SPRAGUE RAPPAPORT	UPC NUMBER
BLACK (01)	SH-101-LS	0-93764-10101-4	DH-201-LS	0-93764-20201-8	SR-301-LS	0-93764-30301-2
NAVY (05)	SH-105-LS	0-93764-10105-2	DH-205-LS	0-93764-20205-6	SR-305-LS	0-93764-30305-0
PURPLE (09)	SH-109-LS	0-93764-60009-8	DH-209-LS	0-93764-60011-1	SR-309-LS	0-93764-30309-8
HUNTER	SH-111-LS	0-93764-60010-4	DH-211-LS	0-93764-60012-8	SR-311-LS	0-93764-30311-1

## Packaging and Shipping Specifications

Selling Unit (Box)	L x W x H (inches)	WT (lb.)	Individual Item Cube (cubic feet)
All Single Heads	4 <sup>7</sup> / <sub>8</sub> x 1 x 10 <sup>1</sup> / <sub>4</sub>	0.25	0.03
All Dual Heads	4 <sup>7</sup> / <sub>8</sub> x 1 x 10 <sup>1</sup> / <sub>4</sub>	0.30	0.03
All Sprague Rappaports	6 <sup>1</sup> / <sub>2</sub> x 1 <sup>5</sup> / <sub>8</sub> x 10 <sup>1</sup> / <sub>4</sub>	0.70	0.06

Unit Case	L x W x H (inches)	WT (lb.)	Case Cube (cubic feet)
All Single Heads (50 units)	17 x 10 <sup>1</sup> / <sub>2</sub> x 16	17.20	1.65
All Dual Heads (50 units)	17 x 10 <sup>1</sup> / <sub>2</sub> x 16	19.50	1.65
All Sprague Rappaports (30 Units)	17 x 11 x 16	24.10	1.73



1756 Automation Parkway  
San Jose, CA 95131  
www.LifeSourceOnline.com



LifeSource is a trademark of A&D Medical. ©2008 A&D Medical. All rights reserved.



Ask your A&D Sales Representative about other *LifeSource* products.

- Blood Pressure Monitors
- Ambulatory Blood Pressure Monitors
- Nebulizers
- Heart Rate Monitors
- Digital Thermometers
- Personal Scales
- Stethoscopes

For questions about *LifeSource* monitors, call our toll-free *LifeSource* HealthLine at 1-888-726-9966 or visit our web site at [www.LifeSourceOnline.com](http://www.LifeSourceOnline.com)

LIFE:STETHOSCOPE 11/02